



RamSoft

PowerServer

IHE Integration Statement



IHE Integration Statement			
Vendor	Product Name	Version	Release Date
RamSoft Inc	PowerServer, Gateway	4.3	April 19, 2007
This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display	None	
Consistent Presentation of Images	Image Display, Print Composer	None	
Access to Radiology Information	Image Display, Report Reader, Report Repository, External Report Repository Access	None	
Key Image Note	Image Display	None	
Simple Image and Numeric Report	Report Creator, Report Reader, Report Repository, External Report Repository Access, Enterprise Report Repository	None	
Post-Processing Workflow	Image Display	None	
Reporting Workflow	Image Display	None	
Portable Data for Imaging	Portable Media Creator, Image Display, Report Reader, Print Composer, Portable Media Importer	None	
NM Image	Image Display	None	
<u>Internet address for vendor's IHE information:</u>			
http://www.ramsoft.biz/statements/IHE.pdf			
Links to Standards Conformance Statements for the Implementation			

RamSoft

HL7	http://www.ramsoft.biz/statements/HL7%20Conformance.pdf	
DICOM	http://www.ramsoft.biz/statements/DICOM%20Conformance.pdf	
Links to general information on IHE		
In Europe: http://www.rsna.org/ihe	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j